

# ***Total Carbon Analysis from Intact Soil Cores with Laser-Induced Breakdown Spectroscopy***

**Michael H. Ebinger  
David D. Breshears  
David A. Cremers  
Clifton W. Meyer  
Kelly L. Goddard  
Pat J. Unkefer**

***Los Alamos National Laboratory***

***NETL Carbon Conference  
May 6-8, 2003***

# ***Measurements for Carbon Inventories***

---

**Are there new methods for accurate, rapid carbon measurements?**

- **Carbon concentrations with accuracy and precision equal to or better than current methods**
- **Data in minutes rather than days, weeks, or months**
- **Measurements while in the field**
- **Different scales: plots to landscapes**
- **Adequate sampling at reasonable cost**

# ***Laser-Induced Breakdown Spectroscopy (LIBS)***

---

- **Field portable**
- **Works with soil cores and discrete samples**
- **Rapid total carbon analysis (minutes)**
- **Inexpensive to operate**



# ***Carbon Distribution at Different Scales***

---


- **Sample plots or landscapes?**
- **How many samples are representative?**
- **How much sampling is affordable?**





# ***Spatial Variation of Carbon in Soil***

- **Prairie soil core, 40 cm depth**
- **Reclaimed soil core 48 cm depth (similar soil type)**
- **Concentration data at 5 cm intervals**
- **Dry combustion analysis**



# Soil Heterogeneity



- Scale model of a soil at 1000X mag.
- Carbon distributed discretely, but throughout the matrix
- Homogeneous?

# ***One Solution: LIBS Analysis of Intact Core***

- **Analysis of soil core 1mm at a time**
- **C concentration data at 1mm to 1 cm resolution**
- **Entire core analyzed in ~30 minutes**






# ***LIBS: Schematic Overview***


- **Nd:YAG Laser at 1086 nm**
- **Pulsed at 10Hz**
- **1mm spot size**




**1mm Laser Spark on Soil**



**Eschelle  
Spectrograph  
& Computer**



**Spectrum from Plasma**






# ***LIBS Spectrum (Partial)***




- **Carbon line:**  
**247.89 nm**
- **Increased intensity = increased carbon**
- **Can detect other elements: H to U**




# ***Dry Combustion vs. LIBS***

---




# *An Intact Soil Core*






# ***LIBS Core Measurements***

- **C/Si Ratio for each 1mm or 1 cm step**
- **Plot with depth in profile**
- **Combine LIBS data to specific intervals**




# ***LIBS Data from Intact Soil Core***



# ***LIBS Data from Intact Soil Core***


**Each LIBS point is 12 to 20 measurements**






# LIBS-Dry Combustion Correlation

- LIBS mean intensity correlates with C/Si ratio
- Deviation at high C: change of soil matrix (Si, bulk density)



# LIBS Benefits



- About 40 minutes to collect LIBS data; 10 minutes to collect and prepare core; all in-field
- About 1 week to collect, prepare and analyze DC samples and data in lab
- Uses *spatial variation* of carbon and includes *soil heterogeneity*

# ***New Methods of Carbon Measurements?***

---

**Does LIBS method provide accurate, rapid carbon measurements?**

- **Carbon concentrations with accuracy and precision equal to or better than current methods**
- **Data in minutes rather than days, weeks, or months**
- **Measurements while in the field**
- **Different scales: plots to landscapes**
- **Adequate sampling at reasonable cost**