

Raman Spectroscopy for the On-Line Analysis of Oxidation States of Oxygen Carrier Particles

Victoria Leichner, John Kirtley, Hergen Eilers

Applied Sciences Laboratory, Institute for Shock Physics,
Washington State University,
PO Box 1495, Spokane, WA 99210


E-mail: eilers@wsu.edu

**Pittsburgh, PA
April 10-12, 2018**

DOE/NETL: FE0027840

ASL/ISP

Chemical Looping Combustion



Goal: Combust fossil fuels in pure O_2 so as to generate pure CO_2 for storage.

Conditions, including:

- **Temperatures: $800\text{ }^\circ\text{C} - 1000\text{ }^\circ\text{C}$**
- **Pressure: $\sim 10\text{ atm}$**
- **Particles constantly moving**

Optimization of process requires ability to identify oxidation state

Oxygen Carrier Particles



<https://www.netl.doe.gov/newsroom/labnotes/labnotes-archive/01-2014>, accessed 9/26/16

This Project:

- Fe₂O₃/Fe₃O₄
- CaSO₄
- CuO

Desired properties include:

- High conversion efficiency
- High reactivity
- Low agglomeration
- Long lifetime
- Low cost
- Low environmental impact

Goal, Objectives, and Vision

Goal:

Develop a sensor for the on-line analysis of the oxidation state of oxygen carrier particles and demonstrate its feasibility.

Objectives:

- (1) Set up and test a Raman spectroscopy system in combination with a pressurized high-temperature sample chamber.**
- (2) Optimize operating parameters of the Raman spectroscopy system and measure the high-temperature spectra of oxygen carriers.**
- (3) Develop an analysis procedure, including statistical modeling and multivariate calibration, for the interpretation of the Raman spectra.**

Long-term Vision:

Monitoring system that can easily be integrated into different types of CLC systems and provide feedback for process control.

Raman Spectroscopy



- **Widely used for the detection/identification of materials.**
- **Demonstrated for standoff/remote single-shot applications.**



<https://www.sciaps.com/raman-spectrometers/>, accessed 9/30/16


Widely used and proven technique.

Raman Spectroscopy



Provides vibrational information unique to material.

Corrections for Raman Spectra




Processing:

- **Instrumental transfer function**
 - Filters
 - Spectrometer
 - Detector
 - Other optical elements
- **Background**
 - Fluorescence
 - Blackbody
 - Cosmic radiation
 - Stray light
 - Laser fluctuation
- **Multi-peak fitting**
 - Peak position
 - FWHM
 - Peak area

Raman spectra require various corrections.


Fitting of Raman Bands



FWHM—thermal broadening

Band center—band migration due to heating

Integrated— calculate temperature from Stokes/antiStokes bands



Peak fitting provides important information for calibration models.

Raman Analysis

- Heat known materials (e.g., Fe_2O_3 , Fe_3O_4) to high temperature (e.g., 800 °C, 900 °C, and 1000 °C) and measure Raman spectra.
- Perform Inverse calibration (determine composition and temperature):

$$x = \alpha_0 + \alpha_1 R_1 + \alpha_2 R_2 + \cdots + \alpha_h R_h$$

$$T = \beta_0 + \beta_1 R_1 + \beta_2 R_2 + \cdots + \beta_k R_k$$

x : Composition (e.g., mol% Fe_2O_3)

T : Temperature


α_i , β_j : fitting parameters

R_j : subsets of the Raman parameters (frequency; FWHM; area)

Yields T and x in the form of linear combinations of the Raman parameters.

1. Li, H., et al. "Feasibility Study of Using High-Temperature Raman Spectroscopy for On-Line Monitoring and Product Control of the Glass Vitrification Process," Energy", PNNL/DOE 1998.
2. Piepel, G. F., et al. "Statistical Modeling of Raman Spectroscopy data from high-temperature glass melts for on-line monitoring of temperature and composition." Quality Engineering 2001, 13, 667-677


Blackbody Radiation at 1000 °C



$$I_{\lambda} = \frac{2hc^2}{\lambda^5 \left(e^{\left(\frac{hc}{\lambda kT} \right)} - 1 \right)}$$

Using 355 nm instead of 785 nm reduces background by more than 5 orders of magnitude.



Blackbody Radiation - Fe_2O_3





Effects of blackbody radiation apparent using 532 nm excitation.

Envisioned Raman Spectroscopy Field Setup

Hot Metal Oxide



Initial Laboratory Setup



**Customized
Heating
Chamber**

Calibration measurements on well-defined samples.

Investigation of OCPs




- **Calcium Sulfate Studies**
 - Pulsed/time gating approach successful for temperatures $>1000^{\circ}\text{C}$
- **Iron Oxide Studies**
 - Pulsed lasers generally not successful because of instability under intense light.
 - CW lasers have proved promising




CaSO₄ – High Temperature Measurements

532 nm



CW, 100 ms




Single shot, long-pulse, 130 μs

Characteristic Raman peaks observed above 1000°C.


CaSO₄ – LIBS

532 nm



Laser induced breakdown spectra (LIBS) observed using laser pulses of sufficient intensity.

$\text{Fe}_2\text{O}_3/\text{Fe}_3\text{O}_4$ – Challenges with Absorption



Ideally, we want high scattering and low absorption.

Fe₂O₃ – Optimizing Light Intensity & Wavelength

Using CW:



	360 nm**	532 nm*	633 nm*	785 nm*
Intensity	$\leq 10^6 \text{ W/cm}^2$	$\leq 10^5 \text{ W/cm}^2$	$\leq 10^5 \text{ W/cm}^2$	$\leq 10^5 \text{ W/cm}^2$
Highest Temperature	1050 °C	700 °C	400 °C	600 °C

**Using hematite powders (212 μm-600 μm)*

***Light intensity only estimated, used densely packed powder*

**Light intensity must be low to avoid LIBS
(creating an advantage for UV excitation).**


Fe₂O₃ – Comparison of Wavelengths



Shorter wavelengths best for avoiding blackbody.


Fe₂O₃ at 1000 °C

Fe₂O₃ spectra at 1000 °C have been successfully collected using 360 nm and 532 nm excitation.



Benchmark high temperature spectrum.


Reference Fe_2O_3 and Fe_3O_4 Spectra at RT




Raman signatures of powders optimized prior to heating mixture sample.

Fe₂O₃/Fe₃O₄ Powder Mixture

Heating



Cooling





Fe₂O₃/Fe₃O₄ can be differentiated up to 600 °C using CW 532 nm.


Next Steps

- **Optimize Collection of Raman Spectra**
 - Further investigate UV Raman
 - Finalize selection of laser wavelength
 - Utilize lock-in amplifier with photomultiplier tube to minimize spectral noise/background
- **Test NETL Samples**
 - Collect reference spectra prior to heating


CuO-Fe₂O₃-Al₂O₃



mag det HV WD spot HFW
250 x ETD 5.00 kV 5.1 mm 3.0 1.19 mm




mag det HV WD spot HFW
1 000 x ETD 5.00 kV 5.1 mm 3.0 298 μm




mag det HV WD spot HFW
5 000 x ETD 5.00 kV 5.1 mm 3.0 59.7 μm

Next Steps

- **Perform Multivariate Statistical Analysis**
 - **Collect reference measurements for calibration**
 - **Test chemometric software for our analysis**
 - **Determine relative mole fraction of OCPs at a given temperature**




Multiple linear regression (MLR) could be a useful chemometric technique for our analysis.



Summary

- **CaSO₄**
 - **Yields good spectra for both cw and low intensity pulses**
 - **Successfully measured spectra above 1000°C**
 - **LIBS observed with high intensity pulses**
- **Fe₂O₃ and Fe₃O₄**
 - **Shorter wavelengths and low intensity light ideal for avoiding LIBS and blackbody radiation**
 - **Benchmark Fe₂O₃ spectrum at 1000 °C achieved**
 - **Fe₂O₃/Fe₃O₄ Raman spectra collected up to 600 °C**
- **Publications/Presentations**
 - **John Kirtley, Victoria Leichner, Benjamin Anderson, Hergen Eilers, “A comparison of pulsed and continuous lasers for high-temperature Raman measurements of anhydrite,” J. Raman Spectrosc. 10.1002/jrs.5356**
 - **John Kirtley, Victoria Leichner, and Hergen Eilers, “Raman spectroscopy of oxygen carrier particles in harsh environments,” Invited Presentation, SPIE-DSS, April 16, 2018, Orlando, FL**

Questions?



Thanks to DOE/NETL: FE0027840

ASL/ISP