

Research as an Experimental Physicist and Optical Engineer at NETL



Gary Lander

NETL Support Contractor



Washington and Jefferson College

Nov. 10, 2023

Disclaimer

This project was funded by the United States Department of Energy, National Energy Technology Laboratory, in part, through a site support contract. Neither the United States Government nor any agency thereof, nor any of their employees, nor the support contractor, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

*Gary R. Lander^{1,2}; Scott Crawford³; Hari Paudel^{3,4}; Jeffrey Wuenschell^{3,4}; Geunsik Lim^{1,2};
Yuhua Duan³; Michael Buric¹*

¹National Energy Technology Laboratory, 3610 Collins Ferry Road, Morgantown, WV 26505, USA

²NETL Support Contractor, 3610 Collins Ferry Road, Morgantown, WV 26505, USA



³National Energy Technology Laboratory, 626 Cochran Mill Road, Pittsburgh, PA 15236, USA

⁴NETL Support Contractor, 626 Cochran Mill Road, Pittsburgh, PA 15236, USA

Optically Detected Magnetic Resonance (ODMR) and Spin Relaxometry Using NV Centers in Nanodiamonds


Motivation

- Quantum sensors are more sensitive than classical sensors
- Nitrogen NV centers in diamond exhibit utilizable quantum properties at room temperature
- Bound electron energy levels and spin relaxation time are dependent on the local environment
- Physical quantities (e.g., magnetic field, electric field, temperature, and strain) affect crystal lattice, which affect energy levels and spin relaxation time
- Measure affects via optically detected magnetic resonance (ODMR) and spin relaxometry
- First, we will use as magnetic field sensor for REE detection





Bound Electron States Depend on Environment

- NV center provides intra-band states
- Splitting between $m_s = +/- 1$ and $m_s = 0$ due to crystal lattice
- Wiggly lines represent photon absorption/emission
- Relaxation from $m_s = +/- 1$ state can take “dark” route
- External magnetic field induces Zeeman splitting of $m_s = +/- 1$ states




Optically-detected Magnetic Resonance

- Continuous wave optical and AC magnetic field excitation
- Electrons excited from $m_s = \pm 1$ have a chance to decay via “dark” route
- When AC B-field resonant with $m_s = 0 \rightarrow m_s = \pm 1$ transition, quenching of fluorescence observed
- Application of external DC B-field induces Zeeman splitting, two resonances can be observed
- Resonance and splitting can be used to measure local B-field




Spin Relaxometry

- Pump-probe-like excitation scheme
- Lone optical pulse acts as pump, polarizes spin
- MW pulse immediately followed by optical pulse acts as probe
- MW pulse modifies spin state, optical pulse probes spin state
- Measure contrast between fluorescence pulses from pump and probe
- Vary time between optical pulses, τ
- Fit contrast vs. τ with exponential, extracted rate is spin relaxation time



ODMR & Spin Relaxometry Schematic





AOM: acousto-optic modulator; M: mirror; L: lens; DBS: dichroic beam splitter; BS: beam splitter; APD: avalanche photodiode; BPF: bandpass filter



Laser-heated Pedestal Growth and Raman Distributed Temperature Sensing System






Motivation

- Optical fiber used in many sensing applications
- Silica fiber cannot withstand harsh environments, single-crystal (SC) fiber can
- SC fiber innately lacks cladding layer, which is necessary to act as efficient optical waveguide
- Developed method to grow SC fiber with cladding layer
- Can be used in Raman distributed temperature sensor (DTS) system in harsh environments







Laser-heated Pedestal Growth

- CO₂ laser source for heating
- “Doughnut” beam shaper surrounds molten zone with light
- Motors advance feedstock (pedestal) and fiber
- Slow process (mm/min)
- Grows pure crystals (no cladding)




Dopant Additions via Regrowth of Sol-gel Coated Fiber


- Grow cladded fibers with two-stage LHPG
 - Sapphire or yttrium aluminum garnet (YAG)
 - Sol-gel (or other) dopant additions
- Evaluate materials' compatibility in energy systems
- Improve fiber performance



Automatic Dopant Segregation through LHPG: Top left: Visible light guiding in GRIN YAG fiber; Top right: EMPA map of Nd concentration in a GRIN YAG fiber; Bottom plots: Co-doped Nd and Ho: YAG fiber dopant concentrations in X (left) and Y (right)





Raman DTS System



- M = mirror
- LCF = laser cleanup filter
- DBS = dichroic beam splitter
- TRF = thermal reject filter
- L = lens
- FUT = fiber under test
- LRF = laser reject filter
- BPF = bandpass filter

- Pulsed ~ 350 ps, 532 nm green laser
- Excites Raman scattering as pulse propagates
- Collects Raman with fast avalanche photodiodes
- Optics designed for sapphire or yttrium aluminum garnet (YAG) fiber
- First interrogator for single-crystal fiber



NETL RESOURCES

VISIT US AT: www.NETL.DOE.gov



@NETL_DOE



@NETL_DOE



@NationalEnergyTechnologyLaboratory

CONTACT:

Gary.Lander@netl.doe.gov

